

MultiPro™ Individual Antibody and Cocktail Staining Protocol

Required Reagents and Consumables

Description	Manufacturer	Part Number
PBS, Ultrapure (1X)	Avantor	K812-500ML
PBS, Ultrapure (10X)	Invitrogen	AM9625
Bovine serum albumin, nuclease-free	Millipore Sigma	126609100GM
Enhanced Blocking Reagent	Proteintech Genomics	G900005
Protector RNase inhibitor	Roche	3335399001
Tween-20, 10%	Thermo Scientific	28320
NP-40, 10%	Thermo Scientific	85124
Paraformaldehyde, 32%	EMS	15714
Human TruStain™ FcX	BioLegend	422302
5mL FACS tubes with caps, nuclease-free	Stemcell Technologies	38057
2mL tubes	Eppendorf	022363352
1.5mL Protein LoBind tubes	Eppendorf	022431081

Required Equipment

- Tabletop centrifuge capable of 14,000 x g, with a rotor suitable for 5 mL FACS tubes.
- Swinging bucket centrifuge with adaptor capable of holding 5mL FACS tubes.
- Tabletop vortex mixer.

MultiPro™ Individual Antibody and Cocktail Staining Protocol

Before You Start

- **Planning your experiment** - Some experiments will require a “Stain-Fix-Stain” approach where live cells are stained with antibodies in between Fc Block and Cell Fixation. If your experiment requires a pre-fixation antibody stain, ensure that you prepare enough BSA Buffer by following the BSA Buffer recipe for experiments requiring a Pre-Fixation Stain.
- **Time required** - Depending on the number of samples being stained, this protocol is expected to take ~3-4 hours *after* cell preparation. The staining protocol *does not* contain any safe stopping points and users must *immediately* proceed to the “Sample Fixation” step of the “[10x Genomics Fixation of Cells & Nuclei for Chromium Fixed RNA Profiling protocol \(CG000478 Rev D\)](#)” and subsequent “Probe Hybridization” of the “[Chromium Fixed RNA Profiling Reagent Kits for Multiplexed Samples with Feature Barcode technology for Protein using Barcode Oligo Capture protocol \(CG000673 Rev A\)](#)”, which requires ~2 hours of hands-on time.
- **Staining with single antibodies** - if you are staining with individual antibodies and *not* using a lyophilized cocktail, skip to the “Sample Preparation” section of this protocol. Make sure to prepare your antibody staining dilution as described in step 23 under the “Intracellular Block & Stain” subsection.
- **Sample preparation** - A high-quality sample is critical for assay performance. The ideal cell suspension has high viability (> 80%) and minimal cell clumps and debris. We strongly recommend straining cells through a 30µm filter such as Pre-Separation Filters (Miltenyi PN: 130-041-407) to remove cell aggregates. During resuspension and mixing steps, pipette gently and avoid introducing bubbles.
- **RNase Awareness** - To prevent RNA degradation, use nuclease-free reagents and consumables listed in this protocol. Before starting, decontaminate your work surfaces and pipettes using RNase Zap (Invitrogen PN: AM9780).
- **Centrifuge selection** – All steps requiring centrifugation of a cell suspension should be performed using a swinging bucket rotor centrifuge. Use of a fixed angle rotor can result in smearing of the cell pellet, which may negatively impact cell recovery.

MultiPro™ Individual Antibody and Cocktail Staining Protocol

Before You Start, *continued*

- **Sequencing Parameters**

For Cell Ranger v8.0 or newer

Read 1: 28 cycles
Index 1: 10 cycles
Index 2: 10 cycles
Read 2: 90 cycles

For Cell Ranger v7.2

Read 1: 48 cycles
Index 1: 10 cycles
Index 2: 10 cycles
Read 2: 50 cycles

MultiPro™ Individual Antibody and Cocktail Staining Protocol

Buffer Preparation

Note: this protocol was developed using human Peripheral Blood Mononuclear Cells (PBMCs). Certain samples may require optimization of detergent and fixative concentrations in the buffers listed below.

Prepare a 30% BSA Stock Solution - Dissolve 30 g of nuclease-free BSA powder in 100 mL of ultrapure water. Filter-sterilize using a 0.2 µm filter and store at 4 °C for up to one week.

BSA Buffer for experiments NOT requiring a Pre-Fixation Stain

Component	[Final]	1 sample (µL)	1 sample +10%	4 samples +10%
RNase inhibitor (200X)	1X	5.0	5.5	22.0
BSA, 30%	1%	33.3	36.6	146.5
10X PBS	1x	100	110	440
Ultrapure water		861.7	947.9	3791.5
Total		1000 µL	1100 µL	4400 µL

BSA Buffer for experiments requiring a Pre-Fixation Stain

Component	[Final]	1 sample (µL)	1 sample +10%	4 samples +10%
RNase inhibitor (200X)	1X	19.5	21.5	85.8
BSA, 30%	1%	130	143	572
10X PBS	1x	390	429	1716
Ultrapure water		3360.5	3696.6	14786.2
Total		3900 µL	4290 µL	17160 µL

Fixation Buffer

Component	[Final]	1 sample (µl)	1 sample + 10%	4 samples +10%
PFA, 32%	4%	25	27.5	110
10X PBS	1X	20	22	88
Ultrapure water	-	155	170.5	682
Total		200 µL	220 µL	880 µL

MultiPro™ Individual Antibody and Cocktail Staining Protocol

IC Staining Buffer

Critical: If you are staining using individual antibodies and not a lyophilized cocktail, use the “single sample + 10%” option when preparing your IC Stain Buffer. Individual antibodies will be diluted together in 55 μ L of IC Staining buffer.

Component	[Final]	1 sample (μ L)	1 sample + 10%	4 samples + 10%
RNase inhibitor (200X)	1x	0.3	0.33	1.3
Tween-20 (10%)	0.2%	1.0	1.1	4.4
NP-40 (10%)	1.0%	5.0	5.5	22.0
BSA (30%)	10%	16.7	18.4	73.5
10X PBS	1x	5.0	5.5	22.0
Ultrapure water	-	22.0	24.2	96.8
Total		50 μL	55 μL	220 μL

IC Blocking Buffer

Component	[Final]	1 sample (μ L)	1 sample + 10%	4 samples + 10%
Enhanced Blocking Reagent	-	15.5	17.1	68.2
RNase inhibitor (200X)	1x	0.3	0.33	1.3
Tween-20 (10%)	0.2%	1.0	1.1	4.4
NP-40 (10%)	1.0%	5.0	5.5	22.0
BSA (30%)	10%	16.7	18.4	73.5
10X PBS	1x	5.0	5.5	22.0
Ultrapure water	-	6.5	7.2	28.6
Total		50 μL	55 μL	220 μL

MultiPro™ Individual Antibody and Cocktail Staining Protocol

IC Wash Buffer

Component	[Final]	1 sample (mL)	1 sample + 10%	4 samples+ 10%
RNase Inhibitor (200X)	1X	0.045	0.050	0.198
Tween-20 (10%)	0.2%	0.180	0.198	0.792
NP-40 (10%)	0.1%	0.090	0.099	0.396
BSA (30%)	1%	0.300	0.330	1.320
10X PBS	1X	0.900	0.990	3.960
Ultrapure Water		7.485	8.234	32.934
Total		9.0 mL	9.9 mL	39.6 mL

MultiPro™ Individual Antibody and Cocktail Staining Protocol

Reconstitution of Lyophilized Antibody Cocktail

1. Remove desired number of cocktail tubes from mylar pouch. If additional tubes remain in the provided pouch, place the extra tube(s) and desiccant pack into a 15 mL conical tube and store at 4 °C.
2. Without removing the cocktail tube cap, equilibrate antibody cocktail tube to room temperature (RT) for 5 minutes.
3. Insert the cocktail tube into a 2 mL Eppendorf tube and spin at 10,000 x g for 30 seconds at RT.

Critical: Step 3 ensures that all lyophilized material is at the bottom of the tube prior to reconstitution. Steps 4 to 7 ensure the reconstituted cocktail is at the correct concentration.

4. Carefully reconstitute the lyophilized cocktail by pipetting the appropriate amount of Staining Buffer down the side of the tube.
 - If using a lyophilized cocktail for pre-fixation staining – 105 µL of BSA Buffer
 - If using a lyophilized cocktail for post-fixation staining – 55 µL of IC Staining Buffer
5. Re-cap the cocktail tube and vortex for 15 seconds and spin at 10,000 x g for 30 seconds at RT.
6. Incubate for 5 minutes at RT.
7. Vortex the tube for 15 seconds and spin at 10,000 x g for 30 seconds at RT.
8. Transfer the entire volume of the reconstituted cocktail into a 1.5 mL Protein LoBind Eppendorf tube and spin at 14,000 x g for 10 minutes at 4 °C.

Critical: Step 8 pellets any protein aggregates that may be present in the reconstituted cocktail.

9. Being careful not to touch the bottom of the tube with your pipette tip, carefully transfer the appropriate volume of reconstituted cocktail to a new 1.5 mL Protein LoBind tube:
 - If using a lyophilized cocktail for pre-fixation staining – 100 µL of BSA Buffer

MultiPro™ Individual Antibody and Cocktail Staining Protocol

- If using a lyophilized cocktail for post-fixation staining – 50 µL of IC Staining Buffer

Critical: Touching the bottom of the tube with a pipette tip or transferring the entire volume could result in the transfer of aggregates to the final staining solution.

10. Keep reconstituted cocktail on ice until ready to use.

MultiPro™ Individual Antibody and Cocktail Staining Protocol

Sample Preparation

1. Prepare 1-2 x 10⁶ cells in 1 mL of BSA Buffer in a 5 ml FACS tube.

Critical: Staining more than 2 x 10⁶ cells per cocktail tube may adversely affect performance.

Fc Block

2. Spin cells at 400 x g for 5 minutes at 4 °C.
3. Using a 1 mL pipette tip, carefully remove supernatant and resuspend cells in 95 µL BSA Buffer. Add 5 µL Human TruStain™ FcX and *gently* pipette mix.

Critical: Centrifugation at low-speed results in loose pellets that can be easily disturbed and aspirated during washing steps. Avoid disruption of the cell pellet as much as possible.

4. Incubate for 10 minutes on ice.
5. Bring up to 1 mL with BSA Buffer.
6. Spin Fc blocked cells at 400 x g for 5 minutes at 4 °C.

Pre-Fixation Stain (optional)

Critical: If omitting Pre-Fixation Staining, skip to step 14.

7. If using single antibodies instead of a lyophilized cocktail - Prepare antibody staining mix by combining antibodies and BSA Buffer to make a total of 100 µL.
8. Remove supernatant from cells and resuspend pellet in 100 µL of BSA Buffer + Antibodies, gently pipette mix.
9. Incubate for 30 minutes on ice.
10. Add 1.4 mL of BSA Buffer, wash cells by gently pipette mixing 10 times.
11. Spin cells at 400 x g for 5 minutes at 4 °C.
12. Remove supernatant

MultiPro™ Individual Antibody and Cocktail Staining Protocol

13. Add 1.4 mL of BSA Buffer, wash cells by gently pipette mixing 10 times.

14. Spin Pre-Fixation-Stained cells at 400 x g for 5 minutes at 4 °C.

Critical: If omitting intracellular staining, remove supernatant and proceed directly to step 28, 10X Genomics sample fixation.

Cell Fixation

15. Being careful not to disturb or aspirate from the cell pellet, remove supernatant.

Critical: Remove as much buffer as possible prior to fixation without disturbing the cell pellet to ensure proper cell fixation, however leaving ~15 µL of residual buffer will not impact assay performance.

16. Set a timer for 15 minutes, then add 200 µL of Fixation Buffer without resuspending the pellet. Immediately start the 15-minute timer and then gently pipette mix *10 times*.

Critical: Carefully follow the sequence of Step 16 to avoid over fixation of cells.

17. Incubate for 15 minutes at room temperature (RT).

18. Bring up to 1 mL with 1X PBS.

19. Spin fixed cell suspension at 850 x g for 5 minutes at 4 °C.

20. Being careful not to disturb or aspirate from the cell pellet, remove supernatant.

Critical: Remove as much Fixation Buffer as possible prior to resuspending cells in IC Block Buffer, however leaving ~15 µL of residual buffer will not impact assay performance.

Intracellular Block & Stain

21. Resuspend cell pellet in 50 µL of IC Block buffer.

22. Incubate on ice for 10 minutes.

MultiPro™ Individual Antibody and Cocktail Staining Protocol

23. If using a lyophilized cocktail, add all 50 µL of reconstituted antibody cocktail and gently pipette mix 10 times.

Critical: If not using a lyophilized cocktail, combine antibody staining volume with enough IC Staining Buffer so that the final volume is 50 µL. Add 50 µL of the diluted antibody pool to the cells in 50 µL of IC Block Buffer.

24. Incubate for 30 minutes on ice.

25. Wash 1: Add 900 µL of IC Wash Buffer to sample. Wash cells by gently pipetting 10 times then add 2 mL of IC Wash Buffer.

26. Spin at 850 x g for 5 minutes at 4 °C.

27. Wash 2: Remove supernatant and resuspend in 1 mL of IC Wash Buffer. Wash cells by gently pipetting 10 times then add 2 mL of IC Wash Buffer.

28. Spin at 850 x g for 5 minutes at 4 °C.

29. Wash 3: Remove supernatant and resuspend in 1 mL of IC Wash Buffer. Wash cells by gently pipetting 10 times then add 2 mL of IC Wash Buffer.

30. Immediately proceed to "Sample Fixation" in the [10x Genomics Fixation of Cells & Nuclei for Chromium Fixed RNA Profiling protocol \(CG000478 Rev D\)](#). After cell fixation using the 10X Genomics protocol, proceed to "Probe Hybridization" in the [Chromium Fixed RNA Profiling Reagent Kits for Multiplexed Samples with Feature Barcode technology for Protein using Barcode Oligo Capture protocol \(CG000673 Rev A\)](#).

31. Post probe hybridization, we recommend the "Individual Wash Workflow" on page 94 of the ["Chromium Fixed RNA Profiling Reagent Kits for Multiplexed Samples with Feature Barcode technology for Protein using Barcode Oligo Capture protocol \(CG000673 Rev A\)"](#) on protocol on day 2 of the protocol. This aids in reducing cell clumping.

Contact Us

C10011 (Rev A) protocol_multipro_staining

MultiPro™ Individual Antibody and Cocktail Staining Protocol

Questions or concerns? Please contact us at genomics.support@ptglab.com