

Flow Cytometry Methanol Permeabilization Protocol

Reagents required:

Methanol, pre-cooled at -20°C

1x PBS

Flow cytometry antibodies

Experiment procedures:

1. Harvest cells and wash them twice with 1x PBS by centrifugation at 350-500 x g for 5 minutes each time, discard the supernatant.
2. Resuspend the cells at a density of 1×10^6 cells in 100 μ L of 90% methanol. Mix well to dissociate the precipitate.
3. Fix and permeate cells for about 15 minutes at room temperature (20-25°C).
4. Centrifuge at 350-500 x g for 5 minutes, discard the supernatant. Wash the cells 3 times with 1x PBS by centrifugation at 350-500 x g for 5 minutes each time.
5. Resuspend the cells at a cell density of approximately 1×10^6 cells in 100 μ L of 3% BSA (or serum), incubate for 30-60 minutes at 4°C in the dark.
6. (Optional) Centrifuge at 350-500 x g for 5 minutes, discard the supernatant. Wash the cells with enough 1x PBS by centrifugation at 350-500 x g for 5 minutes. Resuspend the cells with 1x PBS.
7. Incubate the cells with the primary antibody in each 100 μ L of cell resuspension. The concentration of the primary antibody is based on the recommendations or the results of titration.
8. Incubate for 45-60 minutes at 4°C in the dark.
9. Wash the cells with 1x PBS by centrifugation at 350-500 x g for 5 minutes, discard the supernatant. Repeat.
Note: If using fluorochrome-conjugated primary antibodies, skip to step 13.
10. Resuspend the cells with diluted fluorochrome-conjugated secondary antibody in 100 μ L 1x PBS (use recommended concentration for secondary antibody dilution).
11. Incubate for 45-60 minutes at 4°C in the dark.
12. Wash the cells with 1x PBS by centrifugation at 350-500 x g for 5 minutes, discard the supernatant.
13. Resuspend the cells with 200-500 μ L 1x PBS and analyze on flow cytometer.