

Streptavidin Magnetic Beads

Kit Contents

- 1x streptavidin magnetic beads
(Cat #: **MS001**)



Protocol

1. Take the cells of interest, wash and re-suspend in cell separation buffer - PBS, 0.1% BSA, 2mM EDTA, pH 7.4 (100 μ L for every 10⁷ cells).
2. Incubate cells with optimized concentration of biotin conjugated antibody at 4°C for 30 minutes.
3. Wash cells with 2mL PBS and re-suspend in 100 μ L cell separation buffer.
4. Add streptavidin magnetic beads in each tube and incubate at 4°C for 30 minutes.
5. After incubation add 2mL PBS in the tube and put the tube on magnet rack for 10 minutes.
6. Gently remove supernatant, avoiding contact with the cells bound to magnetic beads.
7. The supernatant contains the depleted cells, the enriched cells remain in the tube.
8. Remove tube from magnet, re-suspend cells in 2mL PBS and wash.
9. Now your cells are ready for further analysis.
10. If required, repeat steps 2-8 on the enriched cells for better results.