For Research Use Only

CoraLite® Plus 488-conjugated P27; KIP1 Polyclonal antibody

Antibodies | ELISA kits | Proteins www.ptglab.com

Catalog Number:CL488-25614


Featured Product

Basic Information	Catalog Number: CL488-25614	GenBank Accession Number: BC001971	Purification Method: Antigen affinity purification					
	Size: 100ul , Concentration: 1000 µg/ml by Nanodrop; Source: Rabbit Isotype: IgG Immunogen Catalog Number: AG22582	GenelD (NCBI): 1027	Recommended Dilutions: IF 1:50-1:500					
		Full Name: cyclin-dependent kinase inhibitor 1B (p27, Kip1) Calculated MW: 198 aa, 22 kDa Observed MW: 27 kDa	Excitation/Emission maxima wavelengths: 493 nm / 522 nm					
				Applications	Tested Applications:	Positive Contr	ontrols:	
					Species Specificity: human, mouse	IF : mouse testis tissue,		
				Background Information	CDKN1B, also named as P27 or KIP1, is a cyclin-dependent kinase inhibitor, which shares a limited similarity with CDK inhibitor CDKN1A/p21. P27 binds to and prevents the activation of cyclin E-CDK2 or cyclin D-CDK4 complexes, and thus controlling cell cycle progression at G1. The degradation of this protein, which is triggered by its CDK dependent phosphorylation and subsequent ubiquitination by SCF complexes, is required for the cellular transition from quiescence to the proliferative state. Downregulation of P27 has been implicated in the progression of several malignancies, including lung cancer, hepatocellular carcinoma, salivary cancer, oral squamous cell carcinomas, and gastric cancer.			
	and thus controlling cell cycle progres dependent phosphorylation and subso from quiescence to the proliferative s malignancies, including lung cancer,	ssion at G1. The degradation of this pro equent ubiquitination by SCF complex :tate. Downregulation of P27 has been	tein, which is triggered by its CDK es, is required for the cellular transition implicated in the progression of several					
Storage	and thus controlling cell cycle progres dependent phosphorylation and subso from quiescence to the proliferative s malignancies, including lung cancer,	ssion at G1. The degradation of this pro equent ubiquitination by SCF complex state. Downregulation of P27 has been hepatocellular carcinoma, salivary car nt. Stable for one year after shipment.	tein, which is triggered by its CDK es, is required for the cellular transition implicated in the progression of several					


For technical support and original validation data for this product please contact:T: 1 (888) 4PTGLAB (1-888-478-4522) (toll freeE: proteintech@ptglab.comin USA), or 1(312) 455-8498 (outside USA)W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed mouse testis tissue using CoraLite® Plus 488 P27; KIP1 antibody (CL488-25614) at dilution of 1:200.



Immunofluorescent analysis of (4% PFA) fixed mouse testis tissue using CoraLite® Plus 488 P27; KIP1 antibody (CL488-25614) at dilution of 1:200.