

For Research Use Only

# HAND1 Polyclonal antibody

Catalog Number: 10882-1-AP



## Basic Information

<b>Catalog Number:</b> 10882-1-AP	<b>GenBank Accession Number:</b> BC021190	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 150ul , Concentration: 133 µg/ml by Bradford method using BSA as the standard;	<b>GeneID (NCBI):</b> 9421	
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> O96004	
<b>Isotype:</b> IgG	<b>Full Name:</b> heart and neural crest derivatives expressed 1	
<b>Immunogen Catalog Number:</b> AG1324	<b>Calculated MW:</b> 24 kDa	

## Applications

**Tested Applications:**  
IF, ELISA

**Species Specificity:**  
human, mouse, rat

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data