

For Research Use Only

ITFG1 Polyclonal antibody

Catalog Number: 11731-1-AP



Basic Information

Catalog Number: 11731-1-AP	GenBank Accession Number: BC024162	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 133 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 81533	Recommended Dilutions: WB 1:1000-1:4000
Source: Rabbit	Full Name: integrin alpha FG-GAP repeat containing 1	
Isotype: IgG	Calculated MW: 612 aa, 68 kDa	
Immunogen Catalog Number: AG2331	Observed MW: 68 kDa	

Applications

Tested Applications: WB,ELISA	Positive Controls: WB : mouse testis tissue,
Species Specificity: human, mouse, rat	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse testis tissue were subjected to SDS PAGE followed by western blot with 11731-1-AP (ITFG1 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.