

For Research Use Only

RECQL Polyclonal antibody

Catalog Number: 14940-1-AP



Basic Information

| | | |
|---|---|--|
| Catalog Number: 14940-1-AP | GenBank Accession Number: BC001052 | Purification Method: Antigen affinity purification |
| Size: 150ul , Concentration: 400 µg/ml by Nanodrop; | GeneID (NCBI): 5965 | Recommended Dilutions: WB 1:1000-1:4000 |
| Source: Rabbit | Full Name: RecQ protein-like (DNA helicase Q1-like) | |
| Isotype: IgG | Calculated MW: 73 kDa | |
| Immunogen Catalog Number: AG6834 | Observed MW: 73 kDa | |

Applications

Tested Applications:
WB, ELISA

Species Specificity:
human, mouse, rat

Positive Controls:

WB : mouse testis tissue, mouse skeletal muscle, mouse kidney

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 14940-1-AP (RECQL antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.