

For Research Use Only

PRKX Polyclonal antibody

Catalog Number: 16651-1-AP



Basic Information

Catalog Number: 16651-1-AP	GenBank Accession Number: BC041073	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 150 µg/ml by Nanodrop and 147 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 5613	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	Full Name: protein kinase, X-linked	
Isotype: IgG	Calculated MW: 358 aa, 41 kDa	
Immunogen Catalog Number: AG10083	Observed MW: 41 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : Jurkat cells, mouse lung tissue
Species Specificity: human, mouse, rat	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Jurkat cells were subjected to SDS PAGE followed by western blot with 16651-1-AP (PRKX antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.