

For Research Use Only

TESK2 Polyclonal antibody

Catalog Number: 17414-1-AP



Basic Information

| | | |
|--|---|--|
| Catalog Number: 17414-1-AP | GenBank Accession Number: BC033085 | Purification Method: Antigen affinity purification |
| Size: 150ul , Concentration: 900 µg/ml by Nanodrop and 400 µg/ml by Bradford method using BSA as the standard; | GeneID (NCBI): 10420 | Recommended Dilutions: WB 1:500-1:1000 |
| Source: Rabbit | Full Name: testis-specific kinase 2 | |
| Isotype: IgG | Calculated MW: 542 aa, 60 kDa | |
| Immunogen Catalog Number: AG10964 | Observed MW: 62 kDa | |

Applications

| | |
|--|---|
| Tested Applications: WB,ELISA | Positive Controls: WB : NCCIT cell, |
| Species Specificity: human, mouse, rat | |

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1%BSA**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



NCCIT cell were subjected to SDS PAGE followed by western blot with 17414-1-AP (TESK2 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.