

For Research Use Only

GRHL1 Polyclonal antibody

Catalog Number: 17644-1-AP



Basic Information

Catalog Number: 17644-1-AP	GenBank Accession Number: BC067519	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 200 µg/ml by Nanodrop and 167 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 29841	Recommended Dilutions: WB 1:200-1:1000
Source: Rabbit	Full Name: grainyhead-like 1 (Drosophila)	
Isotype: IgG	Calculated MW: 618 aa, 70 kDa	
Immunogen Catalog Number: AG11832	Observed MW: 70 kDa	

Applications

Tested Applications: WB,ELISA	Positive Controls: WB : NCCIT cell,
Species Specificity: human, mouse, rat	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



NCCIT cell were subjected to SDS PAGE followed by western blot with 17644-1-AP (GRHL1 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.