

For Research Use Only

BMT2 Polyclonal antibody

Catalog Number: 21744-1-AP



Basic Information

Catalog Number: 21744-1-AP	GenBank Accession Number: BC114615	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 300 µg/ml by Nanodrop;	GeneID (NCBI): 154743	Recommended Dilutions: WB 1:1000-1:6000
Source: Rabbit	Full Name: chromosome 7 open reading frame 60	
Isotype: IgG	Calculated MW: 405 aa, 46 kDa	
Immunogen Catalog Number: AG16476	Observed MW: 46-50 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HeLa cells, U-87 MG cells, mouse liver, rat brain
Species Specificity: human, mouse, rat	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 21744-1-AP (C7orf60 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.