

For Research Use Only

Connexin 45 Polyclonal antibody

Catalog Number: 22184-1-AP

Featured Product



Basic Information

Catalog Number: 22184-1-AP	GenBank Accession Number: BC096213	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 550 µg/ml by Nanodrop;	GeneID (NCBI): 10052	Recommended Dilutions: WB 1:1000-1:4000
Source: Rabbit	Full Name: gap junction protein, gamma 1, 45kDa	
Isotype: IgG	Calculated MW: 396 aa, 45 kDa	
Immunogen Catalog Number: AG16742	Observed MW: 45 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse heart tissue,
Species Specificity: human, mouse	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse heart tissue was subjected to SDS PAGE followed by western blot with 22184-1-AP (Connexin 45 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.