

For Research Use Only

SH2B3 Polyclonal antibody

Catalog Number: 22198-1-AP



Basic Information

Catalog Number: 22198-1-AP	GenBank Accession Number: BC136451	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 350 µg/ml by Nanodrop and 300 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 10019	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	Full Name: SH2B adaptor protein 3	
Isotype: IgG	Calculated MW: 575 aa, 63 kDa	
Immunogen Catalog Number: AG17381	Observed MW: 63 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HEK-293T cells,
Species Specificity: human	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HEK-293T cells were subjected to SDS PAGE followed by western blot with 22198-1-AP (SH2B3 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.