

For Research Use Only

FAM46B Polyclonal antibody

Catalog Number: 23149-1-AP

Featured Product

4 Publications



Basic Information

Catalog Number: 23149-1-AP	GenBank Accession Number: BC014160	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 1000 µg/ml by Nanodrop and 647 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 115572	
Source: Rabbit	Full Name: family with sequence similarity 46, member B	
Isotype: IgG	Calculated MW: 424 aa, 47 kDa	
Immunogen Catalog Number: AG19523		

Applications

Tested Applications:
ELISA

Cited Applications:
IF, IHC, WB

Species Specificity:
human

Cited Species:
human, rat

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Patrick A Fletcher	36093576	Glia	IHC
Tao Liang	30532005	Exp Mol Med	WB
Hongyang Sang	35116892	Transl Cancer Res	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data