

For Research Use Only

LAMB4 Polyclonal antibody

Catalog Number: 25217-1-AP



Basic Information

Catalog Number: 25217-1-AP	GenBank Accession Number: BC140804	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 900 µg/ml by Nanodrop and 380 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 22798	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	Full Name: laminin, beta 4	
Isotype: IgG	Calculated MW: 1761 aa, 194 kDa	
Immunogen Catalog Number: AG19401	Observed MW: 193 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse skin tissue,
Species Specificity: human, rat, mouse	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse skin tissue were subjected to SDS PAGE followed by western blot with 25217-1-AP (LAMB4 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.