

For Research Use Only

PLGRKT Polyclonal antibody

Catalog Number: 26101-1-AP



Basic Information

Catalog Number: 26101-1-AP	GenBank Accession Number: BC008212	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 500 µg/ml by Nanodrop and 307 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 55848	Recommended Dilutions: WB 1:1000-1:4000
Source: Rabbit	Full Name: chromosome 9 open reading frame 46	
Isotype: IgG	Observed MW: 17 kDa	
Immunogen Catalog Number: AG23410		

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : A431 cells, SH-SY5Y cells
Species Specificity: human	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 26101-1-AP (PLGRKT antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.