

For Research Use Only

# KRT71 Polyclonal antibody

Catalog Number: 26149-1-AP



## Basic Information

<b>Catalog Number:</b> 26149-1-AP	<b>GenBank Accession Number:</b> BC103917	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 150ul , Concentration: 1000 µg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 112802	<b>Recommended Dilutions:</b> WB 1:500-1:2000
<b>Source:</b> Rabbit	<b>Full Name:</b> keratin 71	
<b>Isotype:</b> IgG	<b>Observed MW:</b> 57 kDa	
<b>Immunogen Catalog Number:</b> AG23536		

## Applications

<b>Tested Applications:</b> WB, ELISA	<b>Positive Controls:</b> WB : mouse skin tissue,
<b>Species Specificity:</b> human, mouse	

## Background Information

KRT71 belongs to a family of type II keratins that are specifically expressed in the inner root sheath of hair follicles. The type II keratins are clustered in a region of chromosome 12q13.

## Storage


**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



mouse skin tissue were subjected to SDS PAGE followed by western blot with 26149-1-AP (KRT71 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.