

For Research Use Only

PATZ1 Polyclonal antibody

Catalog Number: 26225-1-AP **1 Publications**



Basic Information

Catalog Number: 26225-1-AP	GenBank Accession Number: BC021091	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 1000 µg/ml by Nanodrop and 527 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 23598	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	Full Name: POZ (BTB) and AT hook containing zinc finger 1	
Isotype: IgG	Observed MW: 70-74 kDa	
Immunogen Catalog Number: AG23912		

Applications

Tested Applications: WB,ELISA	Positive Controls: WB : mouse spleen tissue, Raji cells
Cited Applications: IHC, WB	
Species Specificity: human, mouse	
Cited Species: human	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Cuihong Zhao	29926841	Med Sci Monit	WB,IHC

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse spleen tissue were subjected to SDS PAGE followed by western blot with 26225-1-AP (PATZ1 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.