

For Research Use Only

MUC7 Polyclonal antibody

Catalog Number: 26544-1-AP



Basic Information

Catalog Number: 26544-1-AP	GenBank Accession Number: BC025688	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 300 µg/ml by Nanodrop;	GeneID (NCBI): 4589	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	Full Name: mucin 7, secreted	
Isotype: IgG	Calculated MW: 39 kDa	
Immunogen Catalog Number: AG24216	Observed MW: 38-40 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : human saliva,
Species Specificity: human	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human saliva were subjected to SDS PAGE followed by western blot with 26544-1-AP (MUC7 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.