

For Research Use Only

# SH3BP4 Polyclonal antibody

Catalog Number: 29623-1-AP



## Basic Information

<b>Catalog Number:</b> 29623-1-AP	<b>GenBank Accession Number:</b> BC057396	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 150ul , Concentration: 600 µg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 23677	<b>Recommended Dilutions:</b> WB 1:1000-1:8000
<b>Source:</b> Rabbit	<b>Full Name:</b> SH3-domain binding protein 4	
<b>Isotype:</b> IgG	<b>Calculated MW:</b> 963 aa, 108 kDa	
<b>Immunogen Catalog Number:</b> AG31221	<b>Observed MW:</b> 110 kDa	

## Applications

<b>Tested Applications:</b> WB, ELISA	<b>Positive Controls:</b> WB : HepG2 cells, MCF-7 cells, mouse lung tissue, rat lung tissue
<b>Species Specificity:</b> Human, Mouse , Rat	

## Background Information

### Storage


**Storage:**  
Store at -20°C.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 29623-1-AP (SH3BP4 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours.