

For Research Use Only

Evi2a Polyclonal antibody

Catalog Number: 51014-2-AP



Basic Information

Catalog Number: 51014-2-AP	GenBank Accession Number: BC038124	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 250 µg/ml by Nanodrop and 250 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 14017	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	Full Name: ecotropic viral integration site 2a	
Isotype: IgG	Calculated MW: 24 kDa	
Immunogen Catalog Number: AG0478	Observed MW: 26 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse liver tissue,
Species Specificity: mouse	

Background Information

Storage


Storage:
Store at -20°C.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse liver tissue were subjected to SDS PAGE followed by western blot with 51014-2-AP (Evi2a antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.