

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Angptl4 Fusion Protein



Basic Information

Catalog Number:
Ag0543

Size:
50 µg

Form:
Available lyophilized

Species:
rat

Expression Source:
e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:
Not tested

Endotoxin Level:
Please contact the lab for more information

Validated Application:
Blocking peptide

Peptide Sequence:
MRCAPTAGAALVLC AATAGLLSAQGRPAQPEPPRFA
SWDEMNL LAHG LLQLGHGLREHVERTRGQLGALERR
MAACGNACQGP KGTDPKDRVPEGQAPETLQSLQTQ
LKAQNSKIQQ LFQKVAQQQR YLSKQNLRIQNLQSQI
DLLTPTHLDNGVDKTSR GKRLPKMAQLI GLTPNATRLH
RPPRDCQELFQEG ERHSGLFQIQPLGSPFFLVNCEMT
SDGGWTVIQRR LNSVDFNQSW EAYKDFGDPQGE
FWLGLEKMHSITGDRGSQLAVQLQDWDGNAKLLQFP
IHLGGEDTAYS LQLEPTANELGATNVPNGLSLPFST
WDQDHLRGLNCAKSLSGWWFGTCSHSNLNGQ
YFHSIPRQRQRK KGI FWKTWKGRYYP LQATTLIQQ
MEATAAS
(1-405 aa encoded by BC078944)

Selected Validation Data

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) E: proteintech@ptglab.com W: ptglab.com
(toll free in USA), or 1(312) 455-8498
(outside USA)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.