

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

PRKCB Fusion Protein



Basic Information

Catalog Number:

Ag0918

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MADPAAGPPPSEGEESTVRFARKGALRQKNVHEVKN
HKFTARFFKQPTFCSHCTDFIWGFGKQGFQCQVCC
FVVHKRCHEVTFSCPGADKGPASDDPRSKHKFKIHT
YSSPTFCDHCGSLLYGLIHQGMKCDTCMMNVHKRCV
MNVPSLCGTDHTERRGRIYIQAHIDRDVLI VLV RDAKN
LVPMDPNGLSDPYV K L K L I P D P K S E S K Q K T K I K C S L N P
E W N E T F R F Q L K E S D K D R R L S V E I W D W D L T S R N D F M G
S L S F G I S E L Q K A S V D G W F K L L S Q E E G E Y F N V P V P P E G
S E A N E E L R Q K F E R A K I S Q G T K V P E E K T T N T V S K F D N N
G N R D R M K L T D F N F L M V L G K G S F G K V M L S E R K G T D E L Y
A V K I L K K D V V I Q D D D V E C T M V E K R V L A L P G K P
(1-400 aa encoded by BC036472)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.