

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# LDB1 Fusion Protein



## Basic Information

Catalog Number:

Ag10082

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MLDRDVGPTPMYPPTYLEPGIGRHTPYGNQTDYRIFEL  
NKRLQNWTEECDNLWWDFAFTTEFFEDDAMLTITFCL  
EDGPKRYTIGRTLIPRYFRSIFEGGATELYVVKHPKEAF  
HSNFVSLDCDQGSMTQHGKPMFTQVCVEGRLYLE  
FMFDDMMRIKTWHFSIRQHRELI PRSILAMHAQDPQ  
MLDQLSKNITRCGLSNSTLNLYRLCVI  
(1-211 aa encoded by BC000482)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: proteintech@ptglab.com  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.