

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

SLC37A3 Fusion Protein



Basic Information

Catalog Number:

Ag11535

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

GLLVSP EEI G L S G I E A E E N F E E D S H R P L I N G G E N E D E Y E
P N Y S I Q D D S S V A Q V K A I S F Y Q A C C L P G V I P Y S L A Y A C
L K L V N Y S F F W L P F Y L S N N F G W K E A E A D K L S I W Y D V G
G I I G G T L Q G F I S D V L Q K R A P V L A L S L L L A V G S L I G Y S R F
F I G G P S N M I S S A I S A D L G R Q E L I Q R S S E A L A T V T G I V D
G S G S I G A A V G Q Y L V S L I R D K L G W M W V F Y F F I L M T S C T
I V F I S P L I V R E I F S L V L R R Q A H I L R E
(225-441aa encoded by BC046567)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.