

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# INPP5F Fusion Protein



## Basic Information

**Catalog Number:**

Ag12323

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

VKQTKSNVNI~~GNLRKLG~~NFTKPEMKV~~NFLKPNL~~KVNL  
WKS~~DSSLET~~MENTGVMDK~~VQAESD~~GDMSSD~~NDSYH~~  
S~~DEFLT~~NSK~~SDEDRQ~~LANSLESV~~GPIDYV~~LPSCGIIAS  
APRLGSR~~SQSL~~SST~~DSSV~~HAPSEITVAHG~~SGLG~~KGQE  
SPLK~~KSPS~~AGDVHILTGF~~AKPMDI~~YCHRFVQDAQ~~NKV~~  
THLSE~~TRSV~~SQ~~QASQ~~ERNQ~~MTNQ~~V~~SDET~~QSESTE~~Q~~  
TPSR~~PSQL~~DV~~SLSA~~TG~~PQFL~~SVEPAH~~SVA~~SQKTPT~~SA~~  
SS~~MLE~~LETGLH~~VTPSP~~ESS~~SSRA~~VSPFAKIR~~SSM~~VQV  
ASITQAGL~~THGIN~~FAVSKV~~QKSP~~PEPEIINQV~~QQNEL~~  
KKMFI~~QCQ~~TRIIQI

(C-term-345aa encoded by BC111493)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E : proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.