

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

HS3ST5 Fusion Protein



Basic Information

Catalog Number:
Ag12532

Size:
50 µg

Form:
Available lyophilized

Species:
human

Expression Source:
e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:
Not tested

Endotoxin Level:
Please contact the lab for more information

Validated Application:
Blocking peptide

Peptide Sequence:
MLFKQQA~~W~~L~~R~~Q~~K~~L~~L~~V~~L~~G~~S~~L~~A~~V~~G~~S~~L~~L~~L~~Y~~L~~V~~A~~R~~V~~G~~S~~L~~D~~R~~L~~
Q~~P~~I~~C~~P~~I~~E~~G~~R~~L~~G~~G~~A~~R~~T~~Q~~A~~E~~F~~P~~L~~R~~A~~L~~Q~~F~~K~~R~~G~~L~~L~~H~~E~~F~~R~~K~~G~~N~~
A~~S~~K~~E~~Q~~V~~R~~L~~H~~D~~L~~V~~Q~~Q~~L~~P~~K~~A~~I~~I~~I~~G~~V~~R~~K~~G~~G~~T~~R~~A~~L~~L~~E~~M~~L~~N~~L~~H~~
P~~A~~V~~V~~K~~A~~S~~Q~~E~~I~~H~~F~~F~~D~~N~~D~~E~~N~~Y~~G~~K~~G~~I~~E~~W~~Y~~R~~K~~K~~M~~P~~F~~S~~Y~~P~~Q~~
Q~~I~~T~~I~~E~~K~~S~~P~~A~~Y~~F~~I~~T~~E~~E~~V~~P~~E~~R~~I~~Y~~K~~M~~N~~S~~S~~I~~K~~L~~L~~I~~V~~R~~E~~P~~T~~T~~R~~A~~I~~
S~~D~~Y~~T~~Q~~V~~L~~E~~G~~K~~E~~R~~K~~N~~I~~T~~Y~~Y~~K~~F~~E~~K~~L~~A~~I~~D~~P~~N~~T~~C~~E~~V~~N~~T~~K~~Y~~K~~A~~
V~~R~~T~~S~~I~~Y~~T~~K~~H~~L~~E~~R~~W~~L~~K~~Y~~F~~P~~I~~E~~Q~~F~~H~~V~~V~~D~~G~~D~~R~~L~~I~~T~~E~~P~~L~~E~~L~~Q~~
L~~V~~E~~K~~F~~L~~N~~L~~P~~P~~R~~I~~S~~Q~~Y~~N~~L~~Y~~F~~N~~A~~T~~R~~G~~F~~Y~~C~~L~~R~~F~~N~~I~~I~~F~~N~~K~~C~~L~~A
G~~S~~K~~G~~R~~I~~H~~P~~E~~V~~D~~P~~S~~V~~I~~T~~K~~L~~R~~K~~F~~F~~H~~P~~F~~N~~Q~~K~~F~~Y~~Q~~I~~T~~G~~R~~T~~L~~N~~
W~~P~~
(1-346aa encoded by BC093911)

Selected Validation Data

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) E: proteintech@ptglab.com W: ptglab.com
(toll free in USA), or 1(312) 455-8498
(outside USA)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.