

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# MMP3 Fusion Protein



## Basic Information

**Catalog Number:**

Ag12648

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

GEDTSMNLVQKYLENYYDLEKDVKQFVRRKDSGPVV  
KKIREMQKFLGLEVTGKLDSDTLEVMRKPRCGVPDVG  
HFRTFPGIPKWRKTHLYRIVNYTPDLPKDAVDSAVEK  
ALKVWEEVTPLTFSRLYEGEADIMISFAVREHGDFYPF  
DGPGNVLAHAYAPGPGINGDAHFDDDEQWTKDTTG  
TNLFLVAAHEIGHSLGLFHSANTEALMYPLYHSLDLDL  
RFRLSQDDINGIQSLYGPPDSETPPLVPEPVPPEPG  
TPANCDPALSFDAVSTLRGEILIFKDRHFWRKSLRKLEP  
ELHLISSFWPSLPSGVDAAYEVTSKDLVFIFKGNQFW  
AIRGNEVRAGYPRGIHTLGFPPPTVRKIDAAISDKEKNK  
TYFFVEDKYWRFDEKRNSMEPGFPKQIAEDFPGIDSKI  
DAVFEFEGFFYFFTGSSQLEFDPNNAKVTHTLKSNSW  
LNC

(26-477 aa encoded by BC074869)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: proteintech@ptglab.com  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.