

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

CKAP5 Fusion Protein



Basic Information

Catalog Number:

Ag12747

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e. coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

NI L S A L L V L L Q D S L L A T A S S P K F S E L V M K L W R M V R L L P
D T I N S I N L D R I L L D I H I F M K V F P K E K L Q C K S E F P I R T L K T
L L H T L C K L K G P K I L D H L T M I D N K N E S E L E A H L C R M M K H
S M D Q T G S K S D K E T E K G A S R I D E K S S K A K V N D F L A E I F K
K I G S K E N T K E G L A E L Y E Y K K Y S D A D I E P F L K N S S Q F F Q
S Y V E R G L R V I E M E R E G K G R I S T S T G I S P Q M E V T C V P T P
T S T V S S I G N T N G E E V G P S V Y L E R L K I L R Q R C G L D N T K Q
D D R P P L T S L L S K P A V P T V A S S T D M L H S K L S Q L R E S R E Q
H Q H S D L D S N Q T H S S G T V T S S S T A N I D D L K K R L E R I K S
S R K

(1683-2032 aa encoded by BC035554)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.