

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# PGLYRP3 Fusion Protein



## Basic Information

Catalog Number:

Ag12783

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PET28a, with N-terminal 6<sup>His</sup>.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

DGRVYEGVGWNIQGLHTQGYNNISLGI A F F G N K I G S S  
P S P A A L S A A E G L I S Y A I Q K G H L S P R Y I Q P L L L K E E T C L D  
P Q H P V M P R K V C P N I I K R S A W E A R E T H C P K M N L P A K Y  
V I I I H T A G T S C T V S T D C Q T V V R N I Q S F H M D T R N F C D I  
G Y H F

(92-244aa encoded by BC069741)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
[ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.