

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

TBC1D22A Fusion Protein



Basic Information

Catalog Number:
Ag12953

Size:
50 µg

Form:
Available lyophilized

Species:
human

Expression Source:
e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:
Not tested

Endotoxin Level:
Please contact the lab for more information

Validated Application:
Blocking peptide

Peptide Sequence:
MASDGARKQFWKRSNSKLPQSIQHVVYGAQHPPFDPL
LHGTLRSTAKMPTTPVKAKRVSFTQEFESNTSDAWD
AGEDDELLAMAAESLNSKVVMETANRVLNRHSQRQ
GRPTLQEGPGLQQKPRPEAEPPSPSGDLRLVKSVMSES
HTSCPAESASDAAPLQRSQSLPHSATVTLGGTSDPST
LSSALSEREASRLDKFKQLLAGPNTDLEELRRLSWSGI
PKPVRPMTWKLKSGYLPANVDRRPA TLQRKQKEYFAF
IEHYYSRNDEVHQDTYRQIHIDIPRMSPEALILQPKV
TEIFERILFIWAI RHPASGYVQGINDLVTPFFVVFICEY
IEAEVDTVDVSCVPAEVLNIEADTYWCMKLLDGI
QDNYTFAQPGIQMKVKMLEELVSRI DEQVHRHLDQH
EVRYLQFAFRWMNLLMREVLPRCTIRLWDTYQSEPD
GFSHFHLVCAAF LVRWRKEILEEKDFQELLLFLQNL
TAHWDDEDISLLAEAYRLKFAFADAPNHYYK
(1-517aa encoded by BC001292)

Selected Validation Data

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) E: proteintech@ptglab.com W: ptglab.com
(toll free in USA), or 1(312) 455-8498
(outside USA)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.