

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

FASTKD3 Fusion Protein



Basic Information

Catalog Number:

Ag12959

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

PLI I K L G K Y V V R H V P H F T N E E L R R V L E A F I Y F G H H D T F F
T K A L E H R V A A V C L T D P E V V C R V M E Y C S R E L I L S K P I L
N A V A E T F V C Q T E K F S P R Q I S A L M E P F G K L N Y L P P N A S
A L F R K L E N V L F T H F N Y F P P K S L L K L L H S C S L N E C H P V N F
L A K I F K P L F L Q R L Q G K E S H L D T L S R A Q L T Q L F L A S V L E C
P F Y K G P K L L P K Y Q V K S F L T P C C S L E T P V D S Q L Y R Y V K I G
L T N L L G A R L Y F A P K V L T P Y C Y T I D V E I K L D E E G F V L P S T
A N E D I H K R I A L C I D G P K R F C S N S K H L L G K E A I K Q R H L Q L
L G Y Q V V Q I P Y H E I G M L K S R R E L V E Y L Q R K L F S Q N T V H
W L Q E

(C-term-350aa encoded by BC001295)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.