

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# HOXA10 Fusion Protein



## Basic Information

Catalog Number:

Ag1315

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MSCSESPAANSFLVDSLISGRGEAGGGGGAGGGG  
GGGYAHGGVYLPAAADLPYGLQSCGLFPTLGGKRNE  
AASPGSGGGGGLGPAHGYGPSIDLWLDAPRSCR  
MEPPDGPPPPQQQPPPPPPPPQAPQATSCSFAQ  
NIKEESSYCLYDSADKCPKVSATAAELAPFRGPPPDG  
CALGTSSGVVPGYFRLSQAYGTAKGYGSGGGGAQ  
QLGAGPFPAQPPGRGFDLPPALASGSADAARKERAL  
DSPPPTLACGSGGGSQGDDEAHASSAAEELSPAPS  
ESSKASPEKDSLNGSKGENAANWLTAKSGRKRCPYT  
KHQTTLEKEFLNMYLTRRERLEISRSVHLTDRQVKIW  
FQNRMMKL

(1-374aa encoded by BC013971)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E : proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.