

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# CASQ2 Fusion Protein



## Basic Information

**Catalog Number:**

Ag13246

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, T-HIS, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MKRTHLFIIVGIYFLSSCRAEEGLNFPTYDGKDRVVSLS  
EKNFKQVLKKYDLLCLYYHEPVSSDKVTQKQFQLKEIV  
LELVAQVLEHKAIGFVMVDAKKEAKLAKKLGFDDEGS  
LYILKGDRTIEFDGEFAADVLEFLLDLIEDPVEIISKLE  
VQAFERIEDYIKLIGFFKSEDSYYKAFEEAAEHFQPYI  
KFFATFDKGVAKKLSLKMNEVDFYEPFMDEPIAIPNKP  
YTEEELVEFVKEHQRPTRLRRLRPEEMFETWEDDLNGIH  
IVAFAEKSDPDGYEFLEILKQVARDNTDNPDSLILWID  
PDDFPLLVAWYEKTFKIDLFRPQIGVVNVTADSVW  
MEIPDDDDLPTAAELEDWIEDVLSGKINTEDDDEDD  
DDNSDEEDNDDSDDDDE

(1-399 aa encoded by BC02288)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: proteintech@ptglab.com  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.