

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# ISL2 Fusion Protein



## Basic Information

Catalog Number:

Ag13321

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MVDII FHYPF LGAMGDH S K K K PGTAMCVGCGSQIHD  
QFILRVSPDLEWHAACLKAECSQYLDETCTCFVRD  
GKTYCKRDYVRLFGIKAKCQVGFSSSDLV MRARDSV  
YHIECFRCSVCSRQLPGDEFSLREHELLCRADHGLLE  
RAAAGSPRSPGPLPGARGLHLPDAGSGRQPALRPHVH  
KQTEKTRVRTVLNEKQLHLRTRCYAANPRPDALMKE  
QLVEMTGLSPRVIRVWFQNKRCDDKKSILMKQLQQ  
QQHSDKTSLQGLTGPLVAGSPIRHENAVQGSAVEV  
QTYQPPWKALSEFALQSDLDQPAFQQLVFSFESGSL  
GNSSGSDVTSLSQLPDTPNMSMVPSPVET  
(1-359aa encoded by BC012136)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: proteintech@ptglab.com  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.