

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# LILRA2 Fusion Protein



## Basic Information

**Catalog Number:**

Ag13366

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e. coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MTPILTVLI CLGSLGPRTHVQAGHLPKPTLWAEPGSV  
IIQGSPVTLRCQGS LQAE EYHLYRENKASWVRRIQE  
PGKNGQFPIPSITWEHAGRYHCQY YSHNHSSEYSDPL  
ELVVTGAYSKPTLSALPSPVVT LGGNVT LQCVSQVAF  
DGFILCKEGEDEHPQRLNSHSHARGWSWAIFSVGPV  
SPSRRWYSYRCYAYDSNSPYVW SLPDLELLVPGVSKK  
PSLSVQPGPMVAPGESLT LQCVSDVGYDRFVLYKEGE  
RDFLQRPGWQPQAGLSQANFT LGPVSPSHGGQYRC  
YSAHNLSSEWSAPSDPLDILITGQFYDRPSLSVQPVPT  
VAPGKNVTLLCQSRGQFHTFL LTKEGAGHPPLHRSE  
HQAQQNQAEFRMGPV TSAHVGT YRCYSSLS  
(N-term-400aa encoded by BC017412)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: proteintech@ptglab.com  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.