

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# CNTROB Fusion Protein



## Basic Information

Catalog Number:

Ag13447

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PET28a, with N-terminal 6\*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

RVELVERLQAMLQAHWDEANQLLSTLPPNPPAPP  
AGPSSPGPQPEPEKEERRVWTMPPMAVALKPVLQQSR  
EARDELPGAPPVLCSSSDLSLLGSPFQSQHSFQPLE  
PKPDLTSSTAGAFSALGAFHPDHRAERPFPEEDPGPD  
GEGLLKQGLPPAQLEGLKNFLHQLLETVPQNNENPSV  
DLLPPKSGPLTVPSWEEAPQVPRIPPVHKTIVPLAM  
ASSLFRVPEPPSSHSQGGSPSSGSPERGGDGLTFPRQ  
LMEVSQLLRLYQARGWGALPAEDLLLYLKRLEHSGTDG  
RGDNVPRRNTDSRLGEIPRKEIPSQAVPRLATAPKTE  
KPPARKKSGHPAPSSMRSGGVWR  
(546-903aa encoded by BC021134)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.