

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# TBC1D10C Fusion Protein



## Basic Information

**Catalog Number:**

Ag13490

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PET28a, with N-terminal 6<sup>+</sup>His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MAQALGEDLVQPPPELQDDSSSLGSDSELSGPGPYRQ  
ADRYGFIGGSSAEPGPGHPPADLIRQREMKWVEMTS  
HWEKTMSSRYKKVKMQCRKGIPSALRARCWPLLCGA  
HVCQKNSPGTYQELAEAPGDPQWMETIGRDLHRQF  
PLHEMFVSPQGHGQQGLLQVLKAYTLRPEQGYCQA  
QGPVAAVLLMHLPPPEAFWCLVQICEVYLPGYGPH  
MEAVRLDAEVFMALLRRLPHVHKHLQQVGVGPLYL  
PEWFCLFARSLPFTVLRVWDAFLSEGARVLFVGLT  
LVRLALGTAEQRGACPGLETLGALRAIPPAQLQEEAF  
MSQVHSVLSERDLQREIKAQLQLPDSAPGPPRP  
QVRLAGAQAIFEAQQLAGVRRGAKPEVPRIVVQPPE  
EPRPPRRKQTRGKTFHGLLTRARGPPIEGPPRPQRGS  
TSFLDTRF

(1-446aa encoded by BC062999)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
[ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.