

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

SERPINA10 Fusion Protein



Basic Information

Catalog Number:

Ag13618

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MKVVPSLLSVLLAQVWLVPLGAPSPQSPETPAPQN
QTSRVVQAPKEEEDEQEASEEKA S EEEKAWLMASR
QQ LAKETSNFGFSLR KISM RHDGNMVFSPFGMSLA
MTGLMLGATGPTETQIKRGLHLQALKPTKPGLLPSLFK
GLRETL SRNLELGLTQGSFAFIHKDFDVKETFNLSKRY
FDTECVPMNFRNASQAKRLMNHYINKETRGIKPKLFD
EINPETKLILVDYILFKGKWLTPFDPVFTEVDTFHLDKY
KTIKVPMMYGAGKFASTFDKNFRCHVLKLPYQGNAT
MLVVLMEKMGDHLALEDYLTDLVETWLRNMKTRNM
EVFFPKFKLDQKYEMHELLRQMGIRRFSPFADLSELS
ATGRNLQVSRVLQRTVIEVDERGTEAVAGILSEITAYS
MPPVIKVD RPFHFMIYEETS GMLLFLGRV VNP TLL
(1-444aa encoded by BC022261)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.