

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

DDRGK1 Fusion Protein



Basic Information

Catalog Number:

Ag13742

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MVAPVWYLVAAALLVGFILFLTRSRGRAASAGQEPLH
NEELAGAGRVAQPGPLEPEEPRA GGRPRRRDLGSRL
QAQRRAQRVAWAEADENEEEA VILAQEEEGVEKPA
ETHLSGKIGAKKLRKLEEKQARKAQREAAEEREERKR
LESQREA EWKKEERLRLEEEQKEEEERKAREEQQR
EH E EY L K L K E A F V V E E E G V G E T M T E E Q S Q S F L T E F I N Y
I K Q S K V L L E D L A S Q V G L R T Q D T I N R I Q D L L A E G T I T G
V I D D R G K F I Y I T P E E L A A V A N F I R Q R G R V S I A E L A Q A S
N S L I A W G R E S P A Q A P A

(1-314 aa encoded by BC000643)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E : proteintech@ptglab.com W:
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.