

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Angptl4 Fusion Protein



Basic Information

Catalog Number:

Ag14117

Size:

50 µg

Form:

Available lyophilized

Species:

rat

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MRCAPTAGAALVLC AATAGLLSAQGRPAQPEPPRFA
SWDEMNL LAHGLLQLGHGLREHVERTRGQLGALERR
MAACGNACQGP KGTDPKDRVPEGQAPETLQSLQTQ
LKAQNSKIQLFQKVAQQQR YLSKQNLRIQNLQSQI
DLLTPTHLDNGVDKTSR GKRLPKMAQLI GLTPNATRLH
RPPRDCQELFQEG ERHSGLFQIQPLGSPFFLVNCEMT
SDGGWTVIQRR LNSVDFNQSW EAYKDFGDPQGE
FWLGLEKMHSITGDRGSQLAVQLQDWDGNAKLLQFP
IHLGGEDTAYS LQLEPTANELGATN VSPNGLSLPFST
WDQDHLRGLNCAKSLSGW WFGTCSHSNLNGQ
YFHSIPRQRQRK KGI FWKTWKGRYYP LQATLLIQP
MEATAAS

(1-405 aa encoded by BC078944)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com W:
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.