

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

SLC37A3 Fusion Protein



Basic Information

Catalog Number:

Ag14305

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

AGGIVIFFGLLVSP E E I G L S G I E A E E N F E E D S H R P L I N G
G E N E D E Y E P N Y S I Q D D S S V A Q V K A I S F Y Q A C C L P G V I
P Y S L A Y A C L K L V N Y
(218-307 aa encoded by BC046567)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E : proteintech@ptglab.com W:
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.