

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# WAS Fusion Protein



## Basic Information

Catalog Number:

Ag1444

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MSGGPMGGRPGGRGAPAVQQNIPSTLLQDHENQRL  
FEMLRKCLTLATAVVQLYLALPPGAEHWTKEHCGA  
VCFVKDNPQKSYFIRLYGLQAGRLLWEQELYSQLVYS  
TPTPFHTFAGDDCQAGLNFADEDEAQAFRALVQEK  
IQKRNQRQSGDRRQLPPPPTPANEERRGGLPPLHP  
GGDQGGPPVGPLSLGLATVDIQNPDISSRYRGLPAP  
GPSPADKKRSGKKKISKADIGAPSGFKHVSHVSWDP  
QNGFDVNNLDPDLRSLFSRAGISEAQLTDAETSKLIYD  
FIEDQGGLEAVRQEMRRQEPPLPPPPSRGGNQLPRP  
PIVGGNKGRSGPLPPVPLGIAPPPPTPRGPPPPGRGGP  
PPPPPPATGRSGPLPPPPPGAGGPPMPPPPPP  
(1-400 aa encoded by BC012738)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.