

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

MED26 Fusion Protein



Basic Information

Catalog Number:
Ag15381

Size:
50 µg

Form:
Available lyophilized

Species:
human

Expression Source:
e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:
Not tested

Endotoxin Level:
Please contact the lab for more information

Validated Application:
Blocking peptide

Peptide Sequence:
QQLDRVDETPGPPHPKPPRCFSFSPRNSRHEGSFAR
QQSLYAPKGSVPSPSPRPQALDATQVPSPLPLAQPST
PPVRRLELLPSAESPVCWLEQPESHQRLAGPGCKAGL
SPAPELLSLAGFSPDSSKADSDAASSGGSDSKKKRYR
PRDYTVNLDGQVAEAGVKPVRLKERKLTDFPMTRQIK
PLTQKEPVRADSPVHMEQQSRTELDKQEAASLQSP
FEQTNWKELSRNEIQSYLSRQSSLLSSGAQTPGAH
HFMSYLLKQEESTRQGARQLHVLVLPQSPPTDLPGLTR
EVTQDDLDRIQASQWPGVNGCQDTQGNWYDWTQ
CISLDPHGDDGRLNILPVYCLD
(251-600 aa encoded by BC110430)

Selected Validation Data

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) E: proteintech@ptglab.com W: ptglab.com
(toll free in USA), or 1(312) 455-8498
(outside USA)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.