

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

TIGD7 Fusion Protein



Basic Information

Catalog Number:

Ag15597

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

DI~~C~~LP~~G~~KK~~I~~K~~I~~N~~E~~R~~L~~S~~A~~F~~L~~C~~A~~N~~A~~D~~G~~T~~H~~K~~L~~S~~I~~I~~I~~G~~K~~S~~K~~L~~P~~K
S~~V~~K~~E~~D~~T~~S~~T~~L~~P~~V~~I~~Y~~K~~P~~S~~K~~D~~V~~W~~F~~T~~R~~E~~L~~F~~S~~E~~W~~F~~F~~Q~~N~~F~~V~~P~~E
V~~R~~H~~F~~Q~~L~~N~~V~~L~~R~~F~~H~~D~~E~~D~~V~~R~~A~~L~~L~~L~~L~~D~~S~~C~~P~~A~~H~~P~~S~~S~~E~~S~~L~~T~~S~~E~~D~~
G~~R~~I~~K~~C~~M~~F~~F~~P~~H~~N~~T~~S~~T~~L~~I~~Q~~P~~M~~N~~Q~~G~~V~~I~~L~~S~~C~~R~~L~~R~~Y~~R~~W~~K~~Q~~L~~E
E~~S~~L~~V~~I~~F~~E~~E~~S~~D~~D~~E~~Q~~E~~K~~G~~D~~K~~G~~V~~S~~K~~I~~K~~I~~Y~~N~~I~~K~~S~~A~~I~~F~~N~~W~~A~~K~~S~~
W~~E~~E~~V~~K~~Q~~I~~T~~I~~A~~N~~A~~W~~E~~N~~L~~L~~Y~~K~~K~~E~~P~~E~~Y~~D~~F~~Q~~G~~L~~E~~H~~G~~D~~Y~~R~~E~~I~~L~~
E~~K~~C~~G~~E~~L~~E~~T~~K~~L~~D~~D~~D~~R~~V~~W~~L~~N~~G~~D~~E~~E~~K~~G~~C~~L~~L~~K~~T~~K~~G~~G~~I~~T~~K~~E~~V
V~~Q~~K~~G~~G~~E~~A~~E~~K~~Q~~T~~A~~E~~F~~K~~S~~A~~V~~R~~E~~S~~L~~D~~Y~~L~~L~~D~~F~~V~~D~~A~~T~~P~~E~~F~~Q~~
R~~F~~H~~F~~T~~L~~K~~E~~M~~Q~~Q~~E~~I~~V~~K~~K~~Q~~F~~Q~~S~~K~~I~~H~~S~~R~~I~~G~~S~~F~~L~~K~~P~~R~~P~~H~~N~~I~~K~~
D~~S~~F~~S~~G~~P~~S~~T~~S~~G~~S~~N~~H

(C-353aa encoded by BC069623)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E : proteintech@ptglab.com W:
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.