

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# RHOT1 Fusion Protein



## Basic Information

**Catalog Number:**

Ag15605

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

DDDLTPEYLFPLKIPPDCTTELNHAYLFQSTFDK  
HDLRDCALSPDELKGLFKVFPYIPWGPVNNNTVCTN  
ERGWITYQGFLSQWLTLYLDVQRCLEYLGYLSILT  
EQESQASAVTVTRDKKIDLQKKQTQRNVFRCNVIGV  
KNC G K S G V L Q A L L G R N L M R Q K K I R E D H K S Y Y A I N T V Y  
V Y G Q E K Y L L L H D I S E S E F L T E A E I I C D V V C L V Y D V S N P K  
S F E Y C A R I F K Q H F M D S R I P C L I V A A K S D L H E V K Q E Y S I  
S P T D F C R K H K M P P P Q A F T C N T A D A P S K D I F V K L T T M A  
M Y P H V T Q A D L K S S T F W L R A S F G A T V F A V L G F A M Y K A  
L L K Q R

(C-342aa encoded by BC060781)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.