

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

C1orf107 Fusion Protein



Basic Information

Catalog Number:

Ag16199

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

LKRPEYEA V F V G N I D D H F R I G V A I L Q R S I R L Y A P F Y S S
D I L I A S P L G L R T I I G G E G E K K R D F D L S S I E L L I I D Q A D I Y L
M Q N W E H V L H L M N H M N L L P L D S H G V D F S R V R M W S L N
N W S K Y Y R Q T L L F G A L Q D A Q I N S V F N K Y C V N M Q G Q V
A V R N V P M T G S I S H V L V Q L P H V F Q R M E A E N L A S V I D A
R F N F F V N K I L P Q Y R D A V M S H T L I Y I P S Y F D F V R L R N Y F K
K E E L N F T H I C E Y T Q K S G V S R A R H F F L Q G E K Q F L L F T E R
F H F Y K R Y T I K G I R N L I F Y E L P T Y P H F Y S E I C N M L R A T N R G
E E A T W T C T V L Y S K Y D A Q R L A A V V G V E R A A Q M L Q S N
K N V H L F I T G E K

(C-350aa encoded by BC022964)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.