

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

RHOT1 Fusion Protein



Basic Information

Catalog Number:

Ag17143

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

DDDLTPEYLFPLKIPPDCTTELNHAYLFQSTFDK
HDLRDCALSPDELKGLFKVFPYIPWGPVNNNTVCTN
ERGWITYQGFLSQWLTLYLDVQRCLEYLGYLSILT
EQESQASAVTVTRDKIDLQKKQTQRNVFRCNVIGV
KNC G K S G V L Q A L L G R N L M R Q K I R E D H K S Y Y A I N T V Y
V Y G Q E K Y L L L H D I S E S E F L T E A E I I C D V V C L V Y D V S N P K
S F E Y C A R I F K Q H F M D S R I P C L I V A A K S D L H E V K Q E Y S I
S P T D F C R K H K M P P P Q A F T C N T A D A P S K D I F V K L T T M A
M Y P H V T Q A D L K S S T F W L R A S F G A T V F A V L G F A M Y K A
L L K Q R

(277-618 aa encoded by BC060781)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.