

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# IRF4 Fusion Protein



## Basic Information

**Catalog Number:**

Ag1762

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e. coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MNLEGGRRGGEFGMSAVSCGNGKLRQWLDQIDSGK  
YPGLVWENEKSI FRIPWKHAGKQDYNREEDAALFKA  
WALFKGKFREGIDKDPPTWKTRLRALNKSNDFEEL  
VERSQLDISDPYKYRIVPEGAKKGAQLTLEDPMMS  
MSHPYTMTPYPSLPAQQVHNYMPPDRSWRDYV  
PDQPHPEIPYQCPMTFGPRGHHWQGPACENGCVT  
GTFYACAPPESQAPGVPTSEPSIRSAEALAFSDCRLHC  
LYYREILVKELTSSPEGCRISHGHTYDASNLDQVLPY  
PEDNGQRKNEKLLSHLARGVVLWMAPDGLYAKRLCQ  
SRIYWDGPLALCNDPRNKLERDQTKLFDTQQFLSEL  
QAFHHGRSLPRFQVTLFCGEEFPDPQRQRKLTAHV  
EPLLARQLYYFAQQNSGHFLRGYDLPEHISNPEDYHRS  
IRHSSIQE

(1-451 aa encoded by BC015752)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.