

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

G6PD Fusion Protein



Basic Information

Catalog Number:

Ag1832

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MAEQVALSRTQVCGILREELFQGD AFHQSDTHIFIM
GASGDLAKKKIYPTIWWLFRDGLLPENTFIVGYARSRL
TVADIRKQSEPFKATPEEKLKLEDFFARNSYVAGQY
DDAASYQRLNSHMNALHLSQANRLFYALPPTVYEA
VTKNIHESCMSQJGWNRIIVEKPFGRDLQSSDRLSNHI
SSLFREDQIYRIDHYLGKEMVQNLMLVLRANRIFGPIW
NRDNIA CVILTFKEPFGTEGRGGYFDEFGIIRDVMQN
HLLQMLCLVAMEKPASTNSDDVRDEKVKVLCISEVQ
ANNVVLGQYVGNPDGEGEATKGYLDDPTVPRGSTTA
TFAAVVLYVENERWDGVPFILLRC GKALNERKAEVRLQ
FHDVAGDIFHQCKRNELVIRVQPNEAVYTKMMTKK
PGMFFNPEESELDTYGNRYKNVKLPDAYERLLDVF
GSMHFVRSDELREAWRIFTPLLHQIELEKPKPIPIYIG
SRGPTEADELMKRVGFQYEGTYKWWNPHKL
(1-515aa encoded by BC000337)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com W:
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.